

Objective	Nanoscribe 3D Large Feature DiLL 2PP 10x/0.3 0000771	Mediginal Medigi	M. M	IND NO 71/XE9 IN THE PROPERTY OF THE PROPERTY	MO PO AT / XED
Mag, WD	10x	20x	25x, 380um	63x, 190um	63x
	resist	resist air lens	resist	resist	resist oil lens
Application	um to mm 3D	2D-2.5D	um to mm 3D,	0.2um to 10's of	New materials,
area			smooth 3D	microns 3D	0.2um to 10's um
NA	0.3	0.5	0.8	1.4	1.4
Micro. Pos.	4	5	1	3	3
Max area sq	825x825	425x425	285x285 sq	140x140	140x140
Max diam	1000	600	400	200	200
Min feature	~2um x,y	~1um x,y	~0.6um x,y	~150nm x,y,	
size (um)	~10um z	~6um z	~3.3um z	~800nm z	
Min spacing between features um	6um	1um	1um	150nm	
Max feature	8mm(IPQ)	Thickness of PR	~3mm (IP-S)	~3mm (IP DIP)	~150um
height (um)			,	,	(w/170um glass)
Typ. Scan	100	2-8	100	100	20
Speed mm/s					
Max write area (mm)		_	25mm ITO: 17mm		_

Objective	Nanoscribe 3D Large Feature DiLL 2Pp 10x/0.3 0000771	UJGON-nelqiqi Olgookokxos rss-aarr	May 1 3 10 8 0 / X 5 8 10 10 10 10 10 10 10 10 10 10 10 10 10	THO HO YL/XE9 GA	Mo no ri/xeg on the state of th
Write speed	6.8 mm <sup>3</sup> /hr, 9min/mm <sup>3</sup>		Shell: .77 mm <sup>3</sup> /hr, 1.3h/mm <sup>3</sup> Solid: .11 mm <sup>3</sup> /hr, 9h/mm <sup>3</sup>	Solid: .0017 mm³/hr 24days/mm³	
Working Dist (um)	700	2100	380	360 Dip In 170 Oil	360 Dip In 170 Oil
Immersion Medium	IP-Q (n=1.487 @ 589nm)	Air	IP-S (n=1.48 @780nm)	IP Dip (n=1.52)	Oil
Resist Resins	IP-Q	AZ, SU8, Ormocorp	IP-S	IP-Dip	IP-L, IP-G
Substrates	Silicon only, does not work with transparent substrates	Silicon, glass	ITO Glass (n=1.624, 700um thick), Si (n=3.71)	Fused Silica (n=1.45), silicon	Borosilicate glass (n=1.52) 170um thick (coverslip)
Slicing (um)	5	3-6 um	1 (0.2-1.0)	0.5 (0.1-0.5)	0.5
Hatch (um)	1	0.7-1.2 um	0.5 (0.2-0.5)	0.2 (0.1-0.5)	0.2
Photoresist	Z=1 (Dill,PR on	Z=1 (Air, PR on	Z=1 (Dill,PR on	Z=1 (Dill,PR on	Z=0 (Oil, PR on
Orientation	bottom)	bottom facing objective)	bottom facing objective)	bottom)	top)





# **Objectives options**



	63x NA 1.4	25x NA 0.8	20x NA 0.5
Structural Characteristics	High 3D resolution	High 3D smoothness	High 2D resolution
Immersion Medium	IP-Dip/Oil	Resin	Air
Working Distance	360 - 170 μm <sup>†</sup>	380 µm	2100 μm
Writing Field	200 μm	400 µm	600 µm
Standard Resin	IP-Dip/IP-L	IP-S	AZ resists
Refractive Index Contrast	>0.04	>0.1	-
Typ. Slicing Distance	500 nm	1000 nm	-
Typ. Hatching Distance	200 nm	300 nm	-

 $<sup>^{\</sup>dagger}\mbox{Working}$  distance depends on configuration and substrate height.

## **Resin and resists options**



## Immersion Configurations

Name	Refractive Index*	Phase	Structure Characteristics	Printing Mode	Configuration
IP-L 780	1.48	Liquid	High resolution	Oil	63x/Cover slip
IP-G 780	1.49	Gel	Flying features	Oil	63x/Cover slip
IP-Dip	1.51	Liquid	High resolution	DiLL	63x/Fused silica
IP-S	1.48	Liquid	High smoothness	DiLL	25x/ITO glass

## Air Configurations

Name	Refractive Index	Phase	Resist Tone
AZ resists	1.6173 - 1.6953 <sup>†</sup>	Solid	Positive
SU-8	1.58	Solid	Negative
Ormocomp	1.52	Liquid	Negative

<sup>\*</sup> Measured at 20°C at 780 nm.

<sup>†</sup> Refractive index depends on the specific resist and wavelength. Unbleached values at 405 nm are shown here.

2PP resin	Refractive index	Phase	Characteristics	Print Set	Objective	Substrate
IP-Dip	1.521	liquid	high resolution	3D SF	63x	fused silica (3D SF DiLL), silicon (3D LF DiLL)
IP-L 780	1.485	liquid	high resolution	3D SF Oil	63x	borosilicate (3D SF Oil)
IP-G 780	1.495 before bake	gel	high resolution, flying features	3D SF Oil	63x	borosilicate (3D SF Oil)
IP-S	1.486	liquid	high smoothness, mesoscale	3D MF	25x	ITO-coated (3D MF DiLL), silicon (3D LF DiLL)
IP-Q	1.487	liquid	mesoscale	3D LF	10x	silicon (3D LF DiLL)
IP-Visio	1.486	liquid	low fluorescence, non- cytotoxic	3D MF	25x	ITO-coated (3D MF DiLL)
AZ 5214E	1.6990	solid	2D lithography	2D ML	20x	silicon (3D LF DiLL)
AZ 9260	1.6963	solid	2D lithography	2D ML	20x	ITO-coated (3D MF DiLL)
AZ MIR 701	1.7039	solid	2D lithography	2D ML	20x	fused silica (3D SF DiLL)
AZ 40XT	1.5851	solid	2D lithography	2D ML	20x	silicon (3D LF DiLL)
SU-8	1.58	solid	2D lithography	2D ML, 3D SF Oil	20x, 63x	fused silica (3D SF DiLL), silicon (3D LF DiLL), borosilicate (3D SF Oil)
Ormocomp	1.520	liquid	transparent	2D ML, 3D SF Oil	20x, 63x	fused silica (3D SF DiLL), silicon (3D LF DiLL), borosilicate (3D SF Oil)

Other 2PP resins used with the Nanoscribe 3D printer.

2PP resin	Refractive index	Phase	Characteristics	Print Set	Objective	Substrate
AZ 1500	1.7123	solid	2D lithography	2D ML	20x	fused silica (3D SF DiLL), silicon (3D LF DiLL)
AZ 6600	1.7112	solid	2D lithography	2D ML	20x	fused silica (3D SF DiLL), silicon (3D LF DiLL)
AZ ECI 3027	1.7014	solid	2D lithography	2D ML	20x	fused silica (3D SF DiLL), silicon (3D LF DiLL)
AZ nLOF	1.6389	solid	2D lithography	2D ML	20x	fused silica (3D SF DiLL), silicon (3D LF DiLL)





Use the 10x immersion objective for printing large features in conjunction with the IP-Q resin.

A resin stop is required (supplied with the 3D LF). This prevents 2PP resin from creeping into the objective nosepiece.

Property	Value	Comments
Objective	10x NA 0.3	see also objective overview
Immersion (DiLL) mode	yes	
Immersion media	2PP resin	IP resins from Nanoscribe are safe to use with the objective
Working distance WD	700 µm	
Objective opening angle a	35°	
Objective lens diameter D	2.5 mm	
Felt ring	not available	use the dedicated resin stop to prevent damage to the printer
Printing field (Galvo Ø)	1000 μm	stage movement allows for a larger effective printing field
Theoretical lateral $(a_{xy})$ and axial resolution $(a_z)$ ; [ref] [aspect ratio]	$a_{xy}=$ 1.6 µm; $a_z=$ 25.4 µm; 16.0	for 780 nm laser wavelength, IP-Q and 20°C; voxel size and line width depend on solution set and print parameters; single lines are very unstable (high aspect ratio; see article for more information)
Standard 2PP resin	IP-Q	
Solution set	3D LF	
Δn required @ 830nm	>0.5	required substrate thickness to distinguish between two interface signals: >700 um (because of the large depth of focus)
Typical slicing distance 5 µm		atomic distriction from ID O
Typical hatching distance	1 μm	standard recipe for IP-Q





The 25x immersion objective is used for printing medium-sized features. The best surface smoothness is achieved in conjunction with IP-S, hence this objective is recommended for printing micro-optics.

A felt ring is available that prevents 2PP resin from creeping into the objective nosepiece.

Two versions of this objective are in use. Each has equivalent optical performance, but differs in shape and felt ring size as well as an additional "SIL" (silicone oil) marking on the adjustment ring for the latest version. The position of the

adjustment ring markers are also mirrored compared to the previous objective revision.

Table 1: Objective properties

Property	Value	Comments
Objective	25x NA0.8	see also objective overview
Immersion (DiLL) mode	yes	
Immersion media	oil/glycol/water/2PP resin	Most IP resins from Nanoscribe are safe to use with this objective
Working distance WD	380 μm	
Objective opening angle a	31°	
Objective lens diameter D	5.2 mm	
Felt ring size	30 mm / 32 mm	depending on objective version
Printing field (Galvo Ø)	400 μm	stage movement allows for a larger printing field
Theoretical lateral $(a_{xy})$ and axial resolution $(a_z);^{ ext{ref}}$ [aspect ratio]	$a_{xy} =$ 595 nm, $a_z =$ 3313 nm; 5.6	for 780 nm laser wavelength, IP-S and 25°C; voxel size and line width depend on solution set and print parameters
Standard 2PP resin	IP-S	
Solution set	3D MF	
Δn required @ 830nm	>0.1	
Typical slicing distance	1 μm	atom day diversion for ID C
Typical hatching distance	0.5 μm	standard recipe for IP-S





The 63x immersion objective is the objective for printing small features. For printing the finest features possible, IP-Dip and IP-L 780 are recommended.

A felt ring prevents 2PP resin or oil creeping into the objective nosepiece.

Property	Value	Comments	
Objective	63x NA1.4	see also objective overview	
Immersion (DiLL) mode	yes		
Immersion medium	oil/2PP resin	with the exception of IP-G (780), IP resins from Nanoscribe are safe to use as immersion medium	
Working distance WD	360 μm (-170 μm)	when working in oil immersion configuration, the range is reduced by 170 $\mu m$ owing to the substrate thickness	
Objective opening angle a	12°		
Objective lens diameter D	5.5 mm		
Felt ring size	22.5 mm		
Printing field (galvo Ø)	200 μm	stage movement allows for a larger printing area	
Theoretical lateral $(a_{xy})$ and axial resolution $(a_z)$ ; [ref] [aspect ratio]	$a_{xy} =$ 340 nm, $a_z =$ 826 nm; 2.4	for 780 nm laser wavelength, IP-Dip and 25°C; voxel size and line width depend on solution set and print parameters	
Standard 2PP resin	IP-Dip	IP-L 780 or IP-G 780 are also suitable	
Solution set	3D SF	or 3D SF Oil	
Δn required @ 830nm	>0.04		
Typical slicing distance	0.3 μm		
Typical hatching distance	0.2 μm	standard recipe for IP-Dip	



The 20x objective is the air objective for 2D and 2.5D lithography. AZ resins give good results in structuring thin- and thick-film resists.

Table 1: Objective properties

Property	Value	Comments	
Objective	20x NA0.5	see also objective overview	
Immersion (DiLL) mode	no	therefore no felt ring required	
Working distance	2100 μm		
Printing field (galvo Ø)	600 μm	stage movement allows for a larger printing field	
Theoretical lateral $(a_{xy})$ and axial resolution $(a_z);^{\lceil \operatorname{ref}  ceil}$ [aspect ratio]	$a_{xy} =$ 951 nm, $a_z =$ 5824 nm; 6.1	for 780 nm laser wavelength and air; voxel size and line width depend on print set and print parameters	
Standard resin	AZ		
Print set	2D ML		
Δn required	no limitation	typically the interface between substrate and air is found	
Typical slicing distance	3-6 µm	standard regine for A7 regins	
Typical hatching distance	0.7-1.2 μm	standard recipe for AZ resins	